

# Download and Register GARDiS Software

**User Manual: UM0115\_3**

## Foreword

Copyright © 2020 TDSi. All rights reserved.

Time and Data Systems International Ltd operate a policy of continuous improvement and reserves the right to change specifications, colours or prices of any of its products without prior notice.

## Guarantee

For terms of guarantee, please contact your supplier.

## Trademarks

Copyright © 2020 Time and Data Systems International Ltd (TDSi).

This document or any software supplied with it may not be used for any purpose other than that for which it is supplied nor shall any part of it be reproduced without the prior written consent of TDSi. Microsoft and Windows are registered trademarks of Microsoft Corporation. All other brands and product names are trademarks or registered trademarks of their respective owners.

## Cautions and Notes

The following symbols are used in this guide:

	<b>CAUTION!</b> This indicates an important operating instruction that should be followed to avoid any potential damage to hardware or property, loss of data, or personal injury.
	<b>NOTE.</b> This indicates important information to help you make the best use of this product.

## Document Control

Issue	Date Issued	Change Summary	Issued by
1	12 <sup>th</sup> September 2018	Initial release	RT
2	14 <sup>th</sup> November 2018	Changed to Product Registration	RT
3	13 <sup>th</sup> October 2020	Updated screenshots and branding	TBA

## Contents

<b>1. Introduction .....</b>	<b>4</b>
<b>2. Registration.....</b>	<b>4</b>
<b>3. Download .....</b>	<b>6</b>
<b>4. Activate Software .....</b>	<b>8</b>
<b>5. Site Set-up.....</b>	<b>9</b>
<b>6. Managing Sites .....</b>	<b>14</b>

## 1. Introduction

This guide will show you how to download and register the GARDiS software. The registration process uses TDSi's **Product Registration** website which, has been developed to improve the product registration process and can be found at <https://tdsi-product-registration.com>.

Key features of the website are:

- Download latest software installers
- Manage products of different sites
- Register serial numbers
- Deactivate products
- Manage and activate extensions to software licensing

## 2. Registration

1. Go to <https://tdsi-product-registration.com/>

Before you are able to login to the TDSi product registration website you will need to register.


To do this click on the **Click here** link.

2. Complete the registration form and click **“Register”**.



The image shows two overlapping form panels. The top panel is titled "Register" and contains the text "Enter your details to register" and a "Registration Form" label. The bottom panel is titled "Login" and contains the text "Enter your details to log in". It features input fields for "Email Address:" and "Password:", a checkbox for "RememberMe", and a link for "Forgotten Password?". A "Login" button is present, along with a link for "Not yet registered? Click here". At the bottom of the login panel, there is a link: "The information you provide will be processed in accordance with TDSi's privacy policy located here".

3. You will be presented with a confirmation window indicating that an email will be sent to the specified email address.



4. Once you have received the email follow the instructions and click on the link. This will return you to the product registration website to complete the registration process.

Complete the form and click on the **“Submit”** button.

## Complete Registration

Please enter the details to complete your registration.

Email Address:	<input type="text"/>
Company Name:	<input type="text"/>
New Password (Leave blank to remain the same):	<input type="text"/>
Confirm Password:	<input type="text"/>
<input type="submit" value="Submit"/>	

You will now be able to log in to the TDSi product registration website.

### 3. Download

Once you've logged into the TDSi product registration site, click the **“Download”** tab.

You will then see the below screen, click the **“GARDiS”** installer link to download it.

You will then see the installer appear in your internet browser's download section.

For example, in Chrome you will see it appear in the bottom left-hand corner. Once downloaded, click the **“installer file”** to begin installing.

Make sure you download the “**GARDiS Installation Manual**” under the **Documents** title also.

### TDSi Product Registration

- Licenses
- Extensions
- Authorise Users
- Users
- Companies
- Sites
- Download**

#### Download

##### GARDiS

[Click here to download latest version of Gardis installer](#)

[Click here to view documentation and Walk-throughs](#)

##### EXgarde

[Click here to download latest version of EXgarde installer](#)

[Click here to download the Track and Trace app](#)

#### Tools

[Click here to download Toolkit](#)

#### Documents

[Gardis installation manual](#)

[EXgarde installation manual](#)

[EXgarde quick start guide](#)

[Toolkit manual](#)

#### Extras

[Microsoft SQL 2014 Management Studio x64](#)

[Microsoft SQL 2014 Management Studio x86](#)

[Microsoft Windows6.1-KB976932-X64](#)

[Microsoft Windows6.1-KB976932-X86](#)

### TDSi Product Registration

- Sites
- Download**

#### Download

##### GARDiS

[Click here to download latest version of Gardis installer](#)

[Click here to view documentation and Walk-throughs](#)

##### EXgarde

[Click here to download latest version of EXgarde installer](#)

[Click here to download the Track and Trace app](#)

#### Tools

[Click here to download Toolkit](#)

#### Documents

[Gardis installation manual](#)

[EXgarde installation manual](#)

[EXgarde quick start guide](#)

[Toolkit manual](#)


#### Extras

[Microsoft SQL 2014 Management Studio x64](#)


[Gardis.zip](#)

## 4. Activate Software

1. Open the **"GARDiS Configuration"** tool.  
C Drive – Program Files (x86) – TDSi – Gardis – Gardis Configuration
2. Right click the **"Gardis.Configuration"** tool and click run as administrator.



3. Click the **"License"** tab then click **"Add License."**



## 5. Site Set-up

1. Click **Sites** link then **Add Site**.



2. Enter the details of the customer site.

**Add Site** [Close]

Site Name:

Address:

Nature Of Business:

Site Contact Name:

Site Contact Email:

- Once your sites have been set up you can add software to each site. To add software to each site, click on the **Registrations** link relevant to the site.

Sites    Download    Language    Help

**Sites**  
[Add Site](#)    Search for serial number    Enter search text    Search

1 to 1 of 1 rows    25    1

Site Name	Address	Contact Name	Contact Email	Company	
<input type="text" value="Enter search text"/>	<a href="#">Registrations</a> <input type="button" value=""/>				
TMS Head Office	1 Percival Road	Tina Oakfield	T.Oakfield@example.com		

1 to 1 of 1 rows    25    1

- In the **Site Registrations** window, click the **“Add Software”** button to begin the process.

## Site Registrations - TMS Head Office, 1 Percival Road

▼  
Software

[Add Software](#)  
You don't have any software activated for this site. To start, click the "Add Software" link above.

▶  
Registered Products

[Add Registered Product](#)  
You don't have any registered products for this site. To start, click the "Add Registered Product" link above.

5. The registration window will appear requesting the Serial Number and the Machine Code. You will find the Serial Number in the email you received after you purchased GARDiS and the Machine Code from the product activation window.

Click **Activate**.



**Activate Software Registration** [X]


**Serial Number:** STD-H7KYUE-6FNJC-2VC3P-RJ32R-GS4QQ-4NQS1

**Machine Code:** TAU0E-2WFU-SWEQ7-NUY7C-MEV7D

**Machine ID:** 1 (Any) ▼

**Activate**

6. A pop up window will appear containing the Activation Code. This will need to be copied into the GARDiS activation window.



**Activation  
Code: STD-  
F81V9-F5M9D-  
R0WGC-CM69E-  
66SQ5-M6E5V**

---

7. Enter your **Name**, **Company** and **Serial Number** into the pop up window then click **Register**.

Gardis

### Register Gardis

You must register to unlock the full capabilities of Gardis, or continue in trial mode to evaluate the software.

**Name** •  
Rob

**Organization**  
TDSi

**Serial Number** •  
STD-32JP78-2G0H8-9SFP2-J33A9-XF894-LNQ

Get Help

Register Try

8. Click **Activate by entering code**.

Gardis

### Activate Gardis

To continue you must activate which verifies that you are using a fully licensed version of the software.

15 of 15 days remaining

Name this machine to help remember where you've installed the software.  
TDSI-WIN7DEV5

→ Activate by Entering a Code

[Try a different serial number](#)


Extend...

Get Help


Some features may be disabled or reduced until you complete activation.  
[More about Activation](#)

Activate Later Cancel

9. Once you have entered your Activation Code from the pop up window, then click the **“Activate”** button.



10. The GARDiS Configuration tool will look like this, once the license has been activated.



## 6. Managing Sites

1. Once you have logged in correctly, you will be presented with the TDSi product registration main page. Here you will be able to view your registered site, search for serial numbers and add additional sites.

The screenshot shows the 'Sites' management interface. At the top, there's a navigation bar with 'Sites', 'Download', 'Language', and 'Help'. Below this, the 'Sites' section includes an 'Add Site' link and a search bar for serial numbers. A table displays the following data:

Site Name	Address	Contact Name	Contact Email	Company
TMS Poole Branch	1 Poole High St	Finlay Beech	F.Beech@example.com	<a href="#">Registrations</a>
TMS Head Office	1 Percival Road	Tina Oakfield	T.Oakfield@example.com	<a href="#">Registrations</a>

2. The registered sites will automatically be displayed when you log in. You can view registered software from each of the sites by clicking on the “**Registrations**” link to the left of each site.

Product Registration | TDSi | x  
 www.tdsi.co.uk/productregistrationeu.html

TDSi Product Registration greg@tdsi.co.uk Log off  
 Sites Language

### Site Registrations - Thomas, Poole

Software

[Add Software](#)

Serial Number	Activation Code		
PRO-NMRTY-4LHT7-JUEKP-EE5FJ-8DJMP-XN5K7	PRO-LQKVD-Q0TVJ-7NPJD-8FWAV-9FWEE-QR1V2	<a href="#">Deactivate</a>	<a href="#">Extensions</a>
EXP-751CXP-GOSEH-324JA-A9PJL-Q0PVJ-7NPJD	EXP-D99MU-1V2U9-H6LPS-STLHY-A9FMV-PLL73	<a href="#">Deactivate</a>	<a href="#">Extensions</a>
ENT-WHSK56-UJWVW8-924JA-A8PJL-2C96AG-TN8X1	ENT-PMDTY-Q2VOD-TPJLQ-69MV9-F8U2X-8CY01	<a href="#">Deactivate</a>	<a href="#">Extensions</a>
ENT-LXUTSL-STLHY-V9FMV-AH3TQ-48H39-WHSK59	ENT-F9MXW-F8AT4-2WATK-2PN2K-LE5FJ-Q0TVJ	<a href="#">Deactivate</a>	<a href="#">Extensions</a>

Registered Products

[Add Registered Product](#)  
 You don't have any registered products for this site. To start, click the "Add Registered Product" link above.

Version: 1.0.4.0 © 2015 - TDSi